

Mirrors can rotate 360 degrees in the pitch direction.

Ideal for applications where the incident light has multiple angles of incidence.

- The angle of the mirror can be fine tuned with the coarse/fine switching clamp.
- The mount is designed to have the reflective surface at the center of rotation of the mount. Mirror thickness does not affect this.
- 0.25 mm fine pitch screw adjusters or differential micrometers can be used to save space and provide finer adjustment for MHAN mounts $\phi 50.8$ mm and under



Guide

- ▶ The RO-20-60 post (diameter $\phi 20$ mm, length 60mm) is included but it can be replaced with other sizes. Special tools are required to remove the post. Different sized post can be specified at the time of purchase.
- ▶ Kinematic mirror holders, MHG-NL, should be used for low optical axes applications. [Reference](#) C014

Attention

- ▶ Beam splitters mounted at 45 degrees will have the beam blocked by the aluminum frame. The BHAN gimbal beamsplitter holders are recommended and have a larger transmitted clear aperture. [Reference](#) C026
- ▶ Use the coarse/fine switching clamp to lock down the mount after the desired adjustment.
- ▶ The post should be well secured before adjusting the mount.

Mirror Mounting Methods

When mounting a mirror in a mirror holder, use gloves or finger cots so that finger prints do not get on the mirror.

When securing a mirror to the gimbal mirror holder, place the reflective surface downward so that the mirror will be tight against the bottom (face side) of the mirror frame. Place a Delrin ring on the mirror from the top, so that it does not scratch the mirror. Secure the retaining ring into the mirror frame using a spanner wrench or similar tool.

Guide:


First, tighten the retaining ring until it just contacts the mirror.

Second, firmly tighten the retaining ring once, until mirror frame and mirror, Delrin ring, and retaining ring are all in tight contact.


Third, loosen the retaining ring until the mirror can move.

Finally, slowly tighten the retaining ring, stopping at the position where the retaining ring is held lightly. So as not to put stress on the mirror.


When shipping or when used in locations with a lot of vibration, it is possible that the retaining ring will come loose, and the mirror will fall off. In this case, either firmly tighten the retaining ring so that it does not come loose, or secure the retaining ring with thread locking adhesive.




How to Use the Coarse / Fine Switching Clamp and Fine Adjustment Screws




Tilt (pitch) coarse movement control



Rotation (Yaw) coarse movement control



Tilt fine movement control




Rotational fine movement control



Outline Drawing

MHAN-S


M6 P1




Part Number	A (mm)	B (mm)	C (mm)	φD (mm)	φE (mm)	G (mm)	F+G (mm)	MAX H (mm)	I (mm)
MHAN-20S	40	56	10	φ17	φ20	30	54	26.5	20.5
MHAN-25.4S	45	66	10	φ22	φ25.4	35	64	27	20.5
MHAN-30S	45	66	10	φ27	φ30	35	64	27	20.5
MHAN-40S	52.5	79.5	12	φ37	φ40	41	76	27.5	20.5
MHAN-50S	60	92	15	φ46	φ50	46	86	29	20.5
MHAN-50.8S	60	92	15	φ47	φ50.8	46	86	29	20.5
MHAN-60S	65	102	15	φ56	φ60	51	96	28.5	20.5

MHAN-30DM

M6 P1




MHA

 MHA-80S: M6 P1
Other: M10 P1


Part Number	A (mm)	B (mm)	C (mm)	φD (mm)	φE (mm)	G (mm)	F+G (mm)	MAX H (mm)	I (mm)
MHA-80S	89	137	20	φ75	φ80	72	133	48	42.5
MHA-100SA	115	177	21	φ95	φ100	101	184	48	45
MHA-130SA	128	205	24	φ124	φ130	116	214	48	45
MHA-150S	140	227	26	φ144	φ150	126	234	48	45

MHAN-50DM

M6 P1



Screw Type

Primary material: Aluminum
Finish: Black Anodized

Part Number	Options specified*	Compatible Optics Diameter [mm]	Compatible Optics Thickness [mm]	Reflected Beam Clear Aperture (45° incidence) [mm]	Fine Adjustment Range Tilt [°]	Fine Adjustment Range Rotation [°]	Fine Adjustment Resolution Tilt [°/rotation]	Fine Adjustment Resolution Rotation [°/rotation]	Weight [kg]
MHAN-20S	N/UU	φ20	2 – 6	φ9.2	±4	±4	about 0.54	about 0.68	0.3
MHAN-25.4S	N/UU	φ25, φ25.4	2 – 6	φ12.7	±4	±4	about 0.54	about 0.68	0.4
MHAN-30S	N/UU	φ30	2 – 6	φ16.3	±4	±4	about 0.54	about 0.68	0.4
MHAN-40S	N/UU	φ40	2 – 8	φ23.3	±4	±4	about 0.45	about 0.55	0.6
MHAN-50S	N/UU	φ50	3 – 11	φ30.4	±4	±4	about 0.35	about 0.48	0.7
MHAN-50.8S	N/UU	φ50.8	3 – 11	φ30.4	±4	±4	about 0.35	about 0.48	0.7
MHAN-60S	N/UU	φ60	3 – 11	φ37.5	±3	±4	about 0.31	about 0.41	0.9
MHA-80S	—	φ80	4 – 15	φ50.9	±3.5	±5	about 0.49	about 0.72	1.6
MHA-100SA	—	φ100	4 – 15	φ65.1	±3.4	±5	about 0.35	about 0.52	1.9
MHA-130SA	—	φ130	7 – 18	φ86.3	±2.9	±4	about 0.30	about 0.42	2.3
MHA-150S	—	φ150	4 – 20	φ100.4	±2.5	±4	about 0.26	about 0.38	2.5

* For specifying options, please refer to "Conversion of Posts, Post Holders and Pedestal Bases of Holders". [Reference](#) C007

Precision Type

Primary material: Aluminum
Finish: Black Anodized

Part Number	Options specified*	Compatible Optics Diameter [mm]	Compatible Optics Thickness [mm]	Fine Adjustment Range Tilt [°]	Fine Adjustment Range Rotation [°]	Fine Adjustment Resolution Tilt [°/rotation]	Fine Adjustment Resolution Rotation [°/rotation]	Ultra Fine Adjustment Resolution Tilt [°/DIV]	Ultra Fine Adjustment Resolution Rotation [°/DIV]	Conversion Tilt [°/DIV]	Conversion Rotation [°/DIV]	Weight [kg]
MHAN-30DM	N/UU	φ30	2 – 6	±4	±4	about 1.08	about 1.35	about 0.11	about 0.14	about 0.002	about 0.002	0.47
MHAN-50DM	N/UU	φ50	3 – 11	±3	±4	about 0.71	about 0.95	about 0.07	about 0.10	about 0.001	about 0.002	0.58

* For specifying options, please refer to "Conversion of Posts, Post Holders and Pedestal Bases of Holders". [Reference](#) C007

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motorized Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Lenses

Prisms

Polarizers

Lasers

Beam Shaping Diffusers

Filters

Shutter

Others

Fiber