

Pellin Broca Prism | PBPQ

RoHS

Catalog
Code

W3410

Application
SystemsOptics &
Optical
CoatingsOpto-
Mechanics

Bases

Manual
StagesActuators &
AdjustersMotorized
StagesLight Sources &
Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

45 Degrees Angle


Retro-reflectors

Equilateral
Dispersing Prisms

Others


Pellin Broca prism is a dispersing Brewster prism and is designed to emit in the direction of perpendicular to the incident beam. When the incident beam from a YAG laser is at Brewster angle it is possible to separate the second harmonic generation beam (532nm) and the fundamental harmonic generation beam (1064nm).

- The Pellin Broca prism is using Brewster angle and the critical angle to reduce reflection losses and obtain high transmittance.
- There is no coating on the surface of the Pellin Broca prism so it can be used with high energy pulsed laser.
- This prism is used at the (Brewster angle) angle of light intensity of the beam of light (invisible) of the YAG fundamental harmonic generation and second harmonic generation beam to minimize reflection by the prism incident surface.
- Make sure to use polarization direction of laser beam parallel to the bottom surface of the prism.
- It can also be used for multi-wavelength oscillation laser spectroscopy of an Argon laser.



Outline Drawing

(in mm)



Chamfer Ridge line
about C0.3
(No chamfer obtuse angle)

Specifications

Material	Synthetic fused silica
Design wavelength	706nm (intermediate of 532nm and 1063nm)
Angle tolerance	<3'
Surface flatness of substrate	$\lambda/10$
Surface Quality (Scratch-Dig)	20-10
Clear aperture	Circle or ellipse inscribed in a rectangular of 90% of the dimension size


Guide

- We can provide Pellin broca prisms custom to the wavelength of your laser upon request.
- Other sizes are available, please contact our Sales Division with your request.


Attention

- Because it deviates from the Brewster angle, the beam emitted from the ultraviolet wavelength is not a non-reflective.
- Although it can also be used as a dispersing prism of non-polarized light, and not allowed to enter in the Brewster angle, it is not emitted at right angles to the incident angle.
- It can also be dispersed incident S polarized laser beam, reflection loss occurs in the incident surface and the exit surface.
- Fingerprints and dirt adhering to the uncoated surface will effect of the total reflection. Please use without touching anything on the uncoated surface.
- A and B dimension is slightly shorter than the actual catalog because it contains chamfer dimension. Dimensional tolerances are defined at the intersection of each side that does not include a chamfer.

Schematic



Equivalent optical system




Angular dispersion of YAG Laser

Brewster wavelength	1064nm	532nm
Incident angle (Brewster angle) [°]	55.39	55.61
Output angle [°]	1064nm	54.93
	532nm	56.30
	355nm	58.09
	266nm	61.01

Typical Transmittance Data

T: Transmission



Specifications

Part Number	A [mm]	B [mm]	C [mm]	θ_1 [°]	θ_2 [°]
PBPQ-30L20-10	30	50	20	56.13	79.50