



Laser line plate mirrors are plate beamsplitters that are optically coated with dielectric multi-layer on the front surface of optical parallels or wedged substrates.

The rear surface is coated with multi-layer anti-reflection.


- Half mirrors divide input beam into reflectance and transmittance ratio of 1:1. A beamsplitter of R:T=1:1 is called "Half Mirror".
- Any loss from the input beams of this product is minimized because dielectric coatings have no absorption properties. However, the input ratio of reflection to transmission depends on wavelength, polarization and angle of incident of input beam.
- Plate beamsplitters have beam deviations on transmission and ghost on rear surface reflections. Wedged substrates are used to prevent ghosting.



Schematic



Optical Parallel



Wedged Substrate

Specifications

Material	BK7, Synthetic fused silica, CaF ₂
Surface Flatness	$\lambda/10$ (PSMH-157 is Polished)
Coating	Front surface: Dielectric multi-layer partial reflection coating Rear surface: Multi-layer anti-reflection coating
Incident angle	45°
Divergence ratio (reflectance : transmittance)	1 : 1
Surface Quality (Scratch-Dig)	10-5 (PSMH-157: 40-20)
Clear aperture	90% of actual aperture

Guide


- ▶ Please contact our Sales Division for customized products. (Customized on size, wavelength or R:T, etc.) [Reference](#) B068
- ▶ We also have ultra-wideband, broadband and cube types.
- ▶ For a guarantee in reflected wavefront error or transmitted wavefront error, please contact our Sales Division with your requests.
- ▶ Wedged types are marked with an arrow on the side of the substrate indicating the thickest point of the wedge.

Attention

- ▶ Should these products not function as a half mirror, please check the polarization characteristics of the light source. Do note that LD laser is linear in polarization.
- ▶ The beam deviation at transmission of a wedged beamsplitter is large compared to a one made of optical parallel.
- ▶ The amount of beam deviation of a beamsplitter depends on the thickness of the substrate and the wavelength or the incident angle of the input beam.
- ▶ Transmission curves are based on actual measurements and may vary with manufacturing lots.
- ▶ Surface flatness is the reflected wavefront distortion of the surface prior to coating.
- ▶ Be sure to wear laser safety goggles when checking optical path and adjusting optical axis.

Outline Drawing

(in mm)



Compatible Optic Mounts

BHAN-30S, -50S / MHG-MP30-NL, MP50-NL

Application
Systems

Optics &
Optical
Coatings

Opto-
Mechanics

Bases

Manual
Stages

Actuators &
Adjusters

MotORIZED
Stages

Light Sources &
Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator

Beam Samplers

Others



Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motoeized Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator

Beam Samplers

Others

Laser Line Plate Half Mirrors | PSMH


Laser Line							
Part Number	Wavelength Range [nm]	Diameter ϕ D [mm]	Thickness t [mm]	Material	Parallelism W	Reflectance:R Transmittance:T (The average value of the P-Polarization and the S-Polarization) [%]	Laser Damage Threshold* [J/cm ²]
PSMH-30C03-P-157	157	ϕ 30	3	CaF ₂	<3'	R=40±10	0.5
PSMH-50C05-P-157	157	ϕ 50	5	CaF ₂	<3'	R=40±10	0.5
PSMH-30C03-10-193	193	ϕ 30	3	Synthetic fused silica	<5"	T=45±5	1
PSMH-30C05-10W-193	193	ϕ 30	5	Synthetic fused silica	1°±5'	T=45±5	1
PSMH-50C05-10-193	193	ϕ 50	5	Synthetic fused silica	<5"	T=45±5	1
PSMH-50C08-10W-193	193	ϕ 50	8	Synthetic fused silica	1°±5'	T=45±5	1
PSMH-30C03-10-248/266	248 – 266	ϕ 30	3	Synthetic fused silica	<5"	T=50±3	2
PSMH-30C05-10W-248/266	248 – 266	ϕ 30	5	Synthetic fused silica	1°±5'	T=50±3	2
PSMH-50C05-10-248/266	248 – 266	ϕ 50	5	Synthetic fused silica	<5"	T=50±3	2
PSMH-50C08-10W-248/266	248 – 266	ϕ 50	8	Synthetic fused silica	1°±5'	T=50±3	2
PSMH-30C03-10-308/355	308 – 355	ϕ 30	3	Synthetic fused silica	<5"	T= Average 50±5	2
PSMH-30C05-10W-308/355	308 – 355	ϕ 30	5	Synthetic fused silica	1°±5'	T= Average 50±5	2
PSMH-50C05-10-308/355	308 – 355	ϕ 50	5	Synthetic fused silica	<5"	T= Average 50±5	2
PSMH-50C08-10W-308/355	308 – 355	ϕ 50	8	Synthetic fused silica	1°±5'	T= Average 50±5	2
PSMH-30C03-10-405	390 – 410	ϕ 30	3	BK7	<5"	T=50±3	2.1
PSMH-30C05-10W-405	390 – 410	ϕ 30	5	BK7	1°±5'	T=50±3	2.1
PSMH-50C05-10-405	390 – 410	ϕ 50	5	BK7	<5"	T=50±3	2.1
PSMH-50C08-10W-405	390 – 410	ϕ 50	8	BK7	1°±5'	T=50±3	2.1
PSMH-30C03-10-1064	1064	ϕ 30	3	BK7	<5"	T=50±3	20
PSMH-30C05-10W-1064	1064	ϕ 30	5	BK7	1°±5'	T=50±3	20
PSMH-50C05-10-1064	1064	ϕ 50	5	BK7	<5"	T=50±3	20
PSMH-50C08-10W-1064	1064	ϕ 50	8	BK7	1°±5'	T=50±3	20

*Laser pulse width 10ns (PSMH-157: 20ns), repetition frequency 20Hz


Typical Reflectance Data & Typical Transmittance Data

R: Reflectance T: Transmission


PSMH-157




PSMH-193




PSMH-248/266




PSMH-308/355



PSMH-405



PSMH-1064



Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

MotORIZED Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator

Beam Samplers

Others