


High Stability Ball Plunger Post Holders

Used for coarse/fine rotation and height adjustments of posts.
The posts has three poits of contact with the post holder as opposed to the conventional two contact points.



- The plunger prevents posts from slipping as well ad smooth adjustment with one hand when the clamp is loosened.
- The adjustment screw on the bottom of the post holder can be rotated in 90 degree increments.
- BRS-12-□□-M6 or BRS-20-□□-M6 series can be directly fixed to a vibration isolator (M6x1).

Guide


- Clamp screws (RS-CB) can be purchased separately.

Attention


- BRS-12-□□-M6 mounting to M16 is not available; for use with M16, use BRS-12-□□ or BRS-20-□□.

Outline Drawing

BRS-12-M6 Series



BRS-20-M6 Series



Inner Diameter $\phi 12$ Primary material: Aluminum Finish: Black Anodized		
Part Number	A [mm]	Weight [kg]
BRS-12-30-M6	30	0.06
BRS-12-40-M6	40	0.07
BRS-12-50-M6	50	0.07
BRS-12-60-M6	60	0.08
BRS-12-80-M6	80	0.1
BRS-12-100-M6	100	0.12
BRS-12-120-M6	120	0.14
BRS-12-130-M6	130	0.15

Inner Diameter $\phi 20$ Primary material: Aluminum Finish: Black Anodized		
Part Number	A [mm]	Weight [kg]
BRS-20-30-M6	30	0.07
BRS-20-40-M6	40	0.09
BRS-20-50-M6	50	0.1
BRS-20-60-M6	60	0.12
BRS-20-80-M6	80	0.16
BRS-20-100-M6	100	0.19
BRS-20-120-M6	120	0.22
BRS-20-130-M6	130	0.24