

OSDVIO-R4-2010M-200t(900H)

Specifications

SUPPORT :

OSDVIO-S-R4-700h

Isolation System

Rigid (Self-Standing type)

Height of Support

700mm

Maximum Load Capacity

1700 kg

Damping

60–70 shore hardness rubber mount

Height Adjustment

+/-20 mm

TABLE TOP :

OSDVIO-B-2010M-200t

Dimensions

2000 x 1000 mm

Thickness

200 mm

Weight

278 Kg

Core

Honeycomb core structure, 0.25 thick foil

Core Cell Size

2.9 cm²

Core Shear Modulus

19339 kgf/cm² (275,000 psi)

Surface Flatness

+/- 0.1 mm over 600 mm²

Top Skin

430 series ferromagnetic stainless steel, 4.0 mm thick

Bottom Skin

4.5 mm thick steel

Side Walls

2.0 mm thick steel plate with highly damped composite wood

Mounting Holes

M6-1.0

Mounting Hole Grid

25 mm grid

Mounting Hole Border

37.5 mm

Mounting Hole Sealing

21 mm deep clean cylindrical cap (chemical resistant nylon material)

Damping

Broadband damping