

RO extension post ROE

RoHS Catalog Code W6115

Stainless steel post with M6 screw at the tip and M6 screw hole at the bottom.
The optical axis can be adjusted by using an extension post (ROE) to the required height without wasting the holder at hand.
It can also be used with post (RO-12-*** or RO-20-***)




Guide

► Custom sizes not listed in the catalog and additional processing. Contact our Sales Division for more information. Use the Question Sheet for Special Order Posts.


Outline Drawing

ROE-12-***



Part Number	A (mm)
ROE-12-20	20
ROE-12-25	25
ROE-12-50	50
ROE-12-100	100
ROE-12-150	150

ROE-20-***



Part Number	A (mm)
ROE-20-20	20
ROE-20-25	25
ROE-20-50	50
ROE-20-100	100
ROE-20-150	150

Specifications				
Primary material: Stainless steel, Finish: None				
Part Number	Price (JP Yen)	Diameter [mm]	Length A [mm]	Weight [kg]
ROE-12-20	1,500	φ12	20	0.017
ROE-12-25	1,500	φ12	25	0.021
ROE-12-50	1,500	φ12	50	0.043
ROE-12-100	1,500	φ12	100	0.088
ROE-12-150	1,500	φ12	150	0.132

Specifications				
Primary material: Stainless steel, Finish: None				
Part Number	Price (JP Yen)	Diameter [mm]	Length A [mm]	Weight [kg]
ROE-20-20	1,900	φ20	20	0.047
ROE-20-25	1,900	φ20	25	0.06
ROE-20-50	1,900	φ20	50	0.121
ROE-20-100	1,900	φ20	100	0.245
ROE-20-150	1,900	φ20	150	0.368