


OSDOTS series (Overhead Table Systems) are good for off-the-isolated table mounting or storage, such as power supplies, controllers, oscilloscopes, and other instruments, so vibration isolation will not be affected. The height of shelf can be easily adjusted and optional hanging under-shelf additionally installed at the desired position under main shelf for more storage or mounting.



- Convenient location of instruments.
- Increased table space.
- mechanical isolation of equipments with fans.
- Overhead shelf adjusted in height, free standing.
- Levelling feet for uneven floors.
- Formed steel welded design, black power coat finished.
- Hanging under-shelf for more instruments (Optional).
- Grounded-electrical outlets (110V / 220V / 60Hz) on each side of the shelf.
- N.F.B (30A) for preventing instruments from overload currency.

Outline Drawing



Specifications

Part Number	Weight [kg]
OSDOTS-4	106
OSDOTS-5	128
OSDOTS-6	135
OSDOTS-7	157
OSDOTS-8	164
OSDOTS-10	193
OSDOTS-12	249

Application Systems

Optics & Optical Coatings

Holders

Bases

Manual Stages

Actuators

MotORIZED Stages

Light Sources

Index

Guide

Vibration Isolation Systems

Optical Tables

Darkrooms/ Dark Boxes

Optical Bases

Posts

Brackets

Adapters

Tools