

Base Clamps M6 Conversion Adapters

BSC
AD-M6

BSC

RoHS

Catalog
Code

W6029


Clamp for use with BSP series M6 base plates.



- Can be used to secure base plates or other flat plates.
- Clamps plate thicknesses of 7 – 10mm.

Outline Drawing

BSC-60/60UU



Guide

- This base clamp allows the user to secure base plates at an arbitrary position about a single tapped hole.

Specifications		Primary material: Aluminum Finish: Black Anodized	
Part Number	A [mm]	B	Weight [kg]
BSC-60	14.5	M6 P1	0.04
BSC-60UU	15.5	1/4-20UNC	0.04

AD-M6

RoHS

Catalog
Code


W6030

M6-M4 Conversion Adapters

- Adapter for converting a mounting hole from M6 to M4.

Outline Drawing

AD-M6-M4




Specifications		Primary material: Aluminum Finish: Black Anodized	
Part Number	Weight [kg]	Quantity [Pieces]	
AD-M6-M4	0.004	10	

M6-M2 Conversion Adapters

- Adapter for converting a mounting hole from M6 to M2.

Outline Drawing

AD-M6-M2



Specifications		Primary material: Aluminum Finish: Black Anodized	
Part Number	Weight [kg]	Quantity [Pieces]	
AD-M6-M2	0.005	10	

Example of Use of Conversion Adapters

