




TrayPXIe

Passive Component Organizer

Coherent Solutions' TrayPXIe stores your fragile passive fiber optic components such as splitters, connector adaptor patchcords, WDM couplers, and isolators in one handy tray in the PXI mainframe.



Key Features

- Built-in fiber tray for easy mounting of components
- 5 or 17 optical ports with connector of your choice
- Protect fragile fiber pigtails of optical components
- Save precious benchtop space
- Organize the optical set-up for higher productivity

Ordering Information

TrayPXIe - XXXX - XX

Configuration number

1001 = 5 connectors (1 slot)

1002 = 17 connectors (1 slot)

Connector Type

FC = FC

SC = SC

Technical Specifications¹

General Specifications	TrayPXIe
Bus Connector	PXIe
PXI slots	1001: 1 slot 1002: 2 slots
Dimensions (HxWxD)	130 mm x 20mm x 215 mm (5.1" x 0.8" x 8.5") 130 mm x 40mm x 215 mm (5.1" x 1.6" x 8.5")
Weight	~ 1 kg (~2.2 lbs)
Operating temperature range	5 °C to 45 °C (41 °F to 113 °F)
Storage temperature range	-40 °C to 70 °C (-40 °F to 158 °F)

To find out more, get in touch with us today.

Coherent Solutions Ltd

Unit A, 28 Canaveral Drive
Rosedale, Auckland 0632
New Zealand

General enquiries: info@coherent-solutions.com
Technical support: support@coherent-solutions.com
Tel: +64 9 478 4849
Fax: +64 9 478 4851

www.coherent-solutions.com



www.linkedin.com/company/coherent-solutions-ltd
www.facebook.com/CoherentSolutionsLtd
www.youtube.com/CoherentSolutionsLtd
www.weibo.com/CoherentSolutionsLtd
i.youku.com/CoherentSolutionsLtd

© 2018 Coherent Solutions Ltd. All rights reserved. No part of this publication may be reproduced, adapted, or translated in any form or by any means without the written permission of Coherent Solutions Ltd.



INTERLAB Sp. z o.o. | Instrumenty badawcze i pomiarowe

Kosiarzy 37 paw. 20, 02-953 Warszawa | +48 22 840 81 80 | interlab@interlab.pl | interlab.pl