

### RO-EE/UU

Catalog Code **W6083**

#### Stainless steel posts.

In case of using this rod, without buying a new rod, it is possible to be set the holder or the like to the various optical axis height.

- The through-hole of the side of the rod is used to insert the hex wrench. By using this, it can be easily attached and detached to the holder or the like.
- Accessories can be selected from the type of M4 and M6 or inch type.



Example


#### Guide

► Custom sizes not listed in the catalog are available. Use the QuestionSheet for Special Order Posts. [WEB Reference](#) Catalog Code **W6053**

#### Outline Drawing

##### RO-EE/UU

special screw SC-M4-14EE...1 screws, CHS-M4-M6-5...1 screws (EE)  
special screw SC-8U-14UU...1 screws, CHS-1/4-8U-5...1 screws (UU)



C: EE type M6P1, UU type 1/4-20UNC  
D: EE type M4, UU type 8-32UNC

#### Specifications

| Part Number     |                 | A<br>[mm] | Weight<br>[kg] |
|-----------------|-----------------|-----------|----------------|
| METRIC          | INCH            |           |                |
| RO-12.7-25.4EE  | RO-12.7-25.4UU  | 25.4      | 0.026          |
| RO-12.7-38.1EE  | RO-12.7-38.1UU  | 38.1      | 0.038          |
| RO-12.7-50.8EE  | RO-12.7-50.8UU  | 50.8      | 0.051          |
| RO-12.7-76.2EE  | RO-12.7-76.2UU  | 76.2      | 0.077          |
| RO-12.7-101.6EE | RO-12.7-101.6UU | 101.6     | 0.102          |
| RO-12.7-152.4EE | RO-12.7-152.4UU | 152.4     | 0.153          |
| RO-12.7-203.2EE | RO-12.7-203.2UU | 203.2     | 0.204          |
| RO-12.7-304.8EE | RO-12.7-304.8UU | 304.8     | 0.306          |

Primary material: Stainless Steel  
Finish: None

### CHS

Catalog Code **W6084**

CHS, as a conversion screw, is processed by different thread at both ends, in case of M6 or 1/4 inch for holder side and M4 or 8-32UNC for stand side are different, it can be used as a connecting screws for mounting.




- We will prepare the two types of screws INCH and METRIC.
- If you have misplaced the included screws for extension rod, it can be purchased separately.

#### Guide

► It is used for extension rod (RO-12.7-\*\*).

#### Outline Drawing

##### CHS



#### Specifications

| Part Number  | A<br>[mm] | B<br>[mm] | C       | D         | Quantity<br>[Pieces] |
|--------------|-----------|-----------|---------|-----------|----------------------|
| CHS-M4-M6-5  | 3         | 5         | M4      | M6 P1     | 10                   |
| CHS-M4-M6-9  | 3         | 9         | M4      | M6 P1     | 10                   |
| CHS-1/4-8U-5 | 1/8 (3.2) | 5         | 8-32UNC | 1/4-20UNC | 10                   |

Primary material: Stainless Steel  
Finish: None

Application  
Systems

Optics &  
Optical  
Coatings

Holders

Bases

Manual  
Stages

Actuators

Motoeized  
Stages

Light  
Sources

Index

Guide

Vibration  
Isolation Systems

Optical Tables

Darkrooms/  
Dark Boxes

Optical Bases

Posts

Brackets

Adapters

Tools